

Mouse EphB4 Recombinant Protein

Catalog Number: PROTP54761

Overview

| | |
|--------------|--|
| Product Name | Mouse EphB4 Recombinant Protein |
| Description | Mouse EphB4 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 16-539aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | Baculovirus |
| Formulation | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol |

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Mouse, P54761, 16-539aa; LEETLLNTKL ETADLKWVTY PQAEGQWEEL SGLDEEQHSV RYEVCDMKR PGGQAHWLRT GWVPRRGAVH VYATIRFTMM ECLSLPRASR SCKETFTVFY YESEADTATA HTPAWMENPY IKVDTVAAEH LTRKRPGAEA TGKVNKTLR LGPLSKAGFY LAFQDQGACM ALLSLHLFYK KCSWLITNLT YFPETVPREL VVPVAGSCVA NAVPTANPSP SLYCREDGQW AEQQVTGCSC APGYEAAESN KVCRCACQGT FKPQIGDESC LPCPANSN NIGSPVCLCR IGYRARS DP RSSPCTTPPS APRSVVHHLN GSTLRLEWSA PLESGGREDL TYAVRCRECR PGGSCPCGG DMFTDPGPRD LVEPVVAIRG LRPDVTYTFE VAALNGVSTL ATGPPPFEPV NVTTDREVPP AVSDIRVTRS SPSSLILSWA IPRAPSGAVL DYEVKYHEKG AEGPSSVRFL KTSENRAELR GLKRGASYLV QVRARSEAGY GPFQGEHHSQ TQLDESESWR EQLAVEHHHH HH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse EphB4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.