

Human recombinant EphB1 protein

Catalog Number: PROTP54762

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant EphB1 protein
Description	Human recombinant EphB1 protein expressed in HEK293 with His-Tag. Sequence domain: 18-540aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P54762, 18-540aa; MEETLMDTRT ATAELGWTAN PASGWEEVSG YDENLNTIRT YQVCNVFEPN QNNWLLTTFI NRRGAHRIYT EMRFTVRDCS SLPNVP GSK ETFNLYYYET DSVIATKKSA FWSEAPYLKV DTIAADESFS QVDFGGRLMK VNTEVRSFGP LTRNGFYLA F QDYGACMSLL SVRVFFKKCP SIVQNFAVFP ETMTGAESTS LVIARGTCIP NAEEDVPIK LYCNGDGEWM VPIGRCTCKP GYEPENSVAC KACPAGTFKA SQEAEGCSHC PSNSRSPAEA SPICTCRTGY YRADFPPEV ACTSVPSGPR NVISIVNETS IILEWHPPRE TGGRDDVTYN

IICKKCRADR RSCSRCDDNV EFVPRQLGLT ECRVSISSLW AHTPYTFDIQ AINGVSSKSP FPPQHVSNI TTNQAAPSTV PIMHQVSATM
RSITLSWPQP EQPNGIILDY EIRYYEKEHN EFNSSMARSQ TNTARIDGLR PGMVYVVQVR ARTVAGYGKF SGKMCFTLT DDDYKSELRE QLP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human recombinant EphB1 protein

For Research Use Only. Not for use in diagnostic procedures.