

## ARSF (NM\_004042) Human Recombinant Protein

Catalog Number: PROTP54793

### Overview

Product Name	ARSF (NM_004042) Human Recombinant Protein
Description	Recombinant protein of human arylsulfatase F (ARSF)
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Concentration

>50 ug/mL as determined by microplate BCA method

### Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

### Amino Acid Sequence

MRPRRPLVFM~~SLV~~CALLN~~TC~~QAHRVHDDKPNIVLIMVDDL~~GIGDL~~G~~CGY~~GN~~DT~~MRTPHIDRLAREGVRLTQHISAASLCSPSRSAFLTGRYP~~IRSGM~~VSSG  
NRRVIQNLAVPAGLPLNETTLAALLKKQGYSTGLIGKWHQGLNCD~~SRSD~~QCHHPYNYGFDYYYGMPFTLV~~DSCWP~~DP~~SRNTE~~LAFESQLWLCVQLVAI  
AILTLTFGKLSG~~WVS~~PWLLIFSMILFIFLLGYAWFSSHTSP~~LYWD~~CLLMRGHEITEQPMKAERAGSIMVKEAISFLERH~~SKET~~FLLFFSFLH~~VHT~~PLPTTDD  
FTGTSKHGLYGD~~NVEEM~~DSMVGKILDAIDDFGLRNNTLVYFTSDHGGHLEARRGHAQLGGWNGIYKGGKGMGGWEGGIRVPGIVRWPGKVPAGRLI  
KEPTSLMDILPTVASVSGGSLPQDRVIDGRDLMP~~LLQGN~~V~~RHSE~~HEFLFHYCGSYLH~~AVRWIP~~KDDSGSVWKAHYVTPVFQPPASGGCYV~~TS~~LCRCFGE  
QVTYHN~~P~~LLFDLSRDPSESTPLTPATEPLYDFVIKKVANALKEHQETIVPVTYQLSELNQGRTWLKPCCGVFPFCLCDKEE~~EV~~SQPRGPNEKR

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



ARSF (NM\_004042) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.