

Human BMP3b/GDF10 Recombinant Protein

Catalog Number: PROTP55107

Overview

Product Name	Human BMP3b/GDF10 Recombinant Protein
Description	Human BMP3b/GDF10 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 369-478aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	10mM Sodium Citrate buffer (pH 3.5) containing 0.1M NaCl, 1mM DTT, 40% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P55107, 369-478aa; MQWDEPRVCS RRYLKVDFAD IGWNEWIISP KSFDAYYCAG ACEFPMPKIV RPSNHATIQS IVRAVGIIPG
IPEPCCVPDK MNSLGVFLFD ENRNWVLKVY PNMSVDTCAC R

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human BMP3b/GDF10 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.