

## Mouse Tyro3/Dtk Recombinant Protein

Catalog Number: PROTP55144

### Overview

Product Name	Mouse Tyro3/Dtk Recombinant Protein
Description	Mouse Tyro3/Dtk Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 31-419aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, P55144, 31-419aa; AGLKLMGAPV KMTVSQGQPV KLNCVSEGMEDPDIHWMKDG TVVQNASQVS ISEHSWIG LLSLKSVERS  
DAGLYWCQVK DGEETKISQS VWLTVGEVVPF FTVEPKDLAV PPNAPFQLSC EAVGPPEPVT IYWWRGLTKV GGPAPSPSVL NVTGVTQRTE  
FSCEARNIKG LATSIPAIVR LQAPPAAPFN TTVTTISSYN ASVAWVPGAD GLALLHSCTV QVAHAPGEWE ALAVVVPVPP FTCLLRNLAP  
ATNYSLRVRC ANALGPSYGDWVVPFQTKGL APARAPQNFH AIRTDSGLIL EWEEVIPEDP GEGPLGPYKL SWVQENGTQD ELMVEGTRAN  
LTDWDPQKDL ILRVCASNAI GDGPWSQPLV VSSH DHAGRQ GPPHSRTSW

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Tyro3/Dtk Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.