

Mouse VE-Cadherin-5 CDH5 Recombinant Protein

Catalog Number: PROTP55284

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse VE-Cadherin-5 CDH5 Recombinant Protein
Description	Mouse VE-Cadherin-5 CDH5 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 25-599aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P55284, 25-599aa; GPNFPQIDTP NMLPAHHRQK RDWIWNQMHI DEEKNESLPH YVGKIKSNVN RQNAKYVLQG EFAGKIFGVD
ANTGNVLAYE RLDREKVSEY FLTALIVDKN TNKNLEQPSS FTVKVHDIND NWPVFSHQVF NASVPEMSAI GTSVIRVTAV DADDPTVAGH
ATVLYQIVKG NEYFSIDNSG LIFTKIKNLD REKQAEYKIV VETQDALGLR GESGTATVMI RLEDINDNFP VFTQSTYTF SVPEDIRVGKP LGFLTVDPPD
EPQNRMTKYS IMQGEYRDTF TIETDPKRNE GIIKPTKSLD YEVIQYTYFY IEATDPTIRY EYLSSTSGKN KAMVTINVL DVEPPVFQRH
FYHFKLPENQ KKPLIGTVVA KDPDKAQRSI GYSIRKTSR QGFFRITKQG NIYNEKELDR ETYAWYNLTV EANELDSRGN PVGKESIVQV

YIEVLDENDN PPEFAQPYEP KVCENAAQGK LVVQISATDK DVVPVNPCKF FALKNEDSNF TLINNHDNTA NITVKYGQFN REHAKFHLYP
VLISDNGVPS LTGTSTLTVG VCKCNEQGEF TFCEEMAAQA GVSIQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Mouse VE-Cadherin-5 CDH5 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.