

## Mouse LRPAP Recombinant Protein

Catalog Number: PROTP55302

### Overview

Product Name	Mouse LRPAP Recombinant Protein
Description	Mouse LRPAP Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 29-360aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P55302, 29-360aa: QPIAGHG GKYSREKNEP EMAAKRESGE EFRMEKLNQL WEKAKRLHLS PVRLAELHSD LKIQERDELN WKKLKVEGLD  
KDGEKEAKLI HNLNVILARY GLDGRKDAQM VHSNALNEDT QDELGDPRLK KLWHKAKTSG KFSSEELDKL WREFLHYKEK IQEYNVLLDT  
LSRAEEGYEN LLSPSDMAHI KSDLTISKHS ELKDRLRSIN QGLDRLRKVS HQGYGSTTEF EEPRVIDLWD LAQSANFTEK ELESFRELK HFEAKIEKHN  
HYQKQLEISH QKLKHVESIG DPEHISRNKE KYVLLLEKTK ELGYKVKKHL QDLSSRVSR RHNEL

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse LRPAP Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.