

Human PKNOX1 Recombinant Protein

Catalog Number: PROTP55347

Overview

Product Name	Human PKNOX1 Recombinant Protein
Description	Human PKNOX1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-436. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 80% by SDS-PAGE

Amino Acid Sequence

Human, P55347, 1-436; MGSSHHHHHH SGLVPRGSH MMATQTLSD SYQDGQMQV VTELKTEQDP NCSEPDAEGV SPPPVESQTP MDVDKQAIYR HPLFPLLALL FEKCEQSTQG SEGTTASFD VDIENFVRKQ EKEGKPFCE DPETDNLNVK AIQVLRHLL ELEKVNELCK DFCSRYIACL KTKMNSETLL SGEPGSPYSP VQSQQIQSAI TGTISPQGIV VPASALQQGN VAMATVAGGT VYQPVTVVTP QGQVVTQTLSPGTIRIQNSQ LQLQLNQDLS ILHQDDGSSK NKRGVLPKHA TNVMRSLWFQ HIGHPYPTED EKKQIAAQTN LTLQVNNWF INARRRILQP MLDSSCETP KTKKKTAQNR PVQRFWPSI ASGVAQPPPS ELTMSGAVV TITTPVNMNV DSLQSLSSDG ATLAVQQVMM AGQSEDESVD STEEDAGALA PAHISGLVLE NSDSLQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PKNOX1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.