

Human recombinant GFR alpha-1 protein

Catalog Number: PROTP56159-1

Overview

Product Name	Human recombinant GFR alpha-1 protein
Description	Human recombinant GFR alpha-1 protein expressed in HEK293 with hlgG-Tag. Sequence domain: 25-424aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5 mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P56159, 25-424aa; DRLDCVKASD QCLKEQSCST KYRTLRCQVA GKETNFSLAS GLEAKDECERS AMEALKQKSL YNCRCKRGMK
KEKNCLRIYW SMYQSLQGND LLEDSPYEPV NSRSLDIFRV VPFISVEHIP KGNNCLDAAK ACNLDDICK YRSAYITPCT TSVSNDVCNR
RKCHKALRQF FDKVPAKHSY GMLFCSCRDI ACTERRRQTI VPVCSYEERE KPNCLNLQDS CKTNYICRSR LADFFTNCQP ESRSVSSCLK
ENYADCLLAY SGLIGTMTP NYIDSSSLSV APWCDCNSG NDLEECLKFL NFFKDNTCLK NAIQAFGNCS DVTVWQPAFP VQTTTATTTT
ALRVKKNKPLG PAGESNEIPT HVLPPCANLQ AQKLSNVSG NTHLCISNGN YEKEGLGASS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant GFR alpha-1 protein

For Research Use Only. Not for use in diagnostic procedures.