

Human WNT3a Recombinant Protein

Catalog Number: PROTP56704

Overview

| | |
|--------------|--|
| Product Name | Human WNT3a Recombinant Protein |
| Description | Human WNT3a Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 19-352aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays. |
| Size | Starting from 20ug |
| Tag | Non-Tagged |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P56704, 19-352aa; MSYPIWWSLA VGPQYSSLGS QPILCASIPG LVPKQLRFCR NYVEIMPSVA EGIKIGIQEC QHQFRGRRWN
CTTVHDSLAI FGPVLDKATR ESAFVHAIAS AGVAFAVTRS CAEGTAAICG CSSRHQGGSPG KGWKWGGCSE DIEFGGMVSR EFADARENRP
DARSAMNRHN NEAGRQAIAS HMHLKCKCHG LSGSCEVKTC WWSQPDFRAI GDFLKDKYDS ASEMVVEKHR ESRGWVETLR PRYTYFKVPT
ERDLVYEEAS PNFCEPNPET GSGFTRDRTC NVSSHGIDGC DLLCCGRGHN ARAERRREKC RCVFHWCCYV SCQECTRVYD VHTCK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human WNT3a Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.