

## Human recombinant BACE-1 protein

Catalog Number: PROTP56817-1

### Overview

Product Name	Human recombinant BACE-1 protein
Description	Human recombinant BACE-1 protein expressed in Baculovirus with His-Tag. Sequence domain: 22-457aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by BCA assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P56817, 22-457aa; TQHGIRLPLR SGLGGAPLGL RLPRETDEEP EEPGRRGSFV EMVDNLRGKS GQGYVEMTV GSPQTLNIL VDTGSSNFAV GAAPHPFLHR YYQRQLSSTY RDLRKGVYVP YTQGKWEDEL GTDLVSIPIHG PNVTVRANIA AITESDKFFI NGSNWEGILG LAYAEIARPD DSLEPFDSL VKQTHVPLNF SLQLCGAGFP LNQSEVLASV GGSMIIGGID HSLYTGSLWY TPIRREWYVE VIVRVEING QDLKMDCKEY NYDKSIVDSG TTNLRLPKKV FEA AVKSIKA ASSTEKFPDG FWLGEQLVCW QAGTTPWNIF PVISLYLMGE VTNQSFRTI LPQQYLRPVE DVATSQDDCY KFAISQSSTG TVMGAVIMEG FYVVFDRARK RIGFAVSACH VHDEFRTAAV EGPFVTL DME DCGYNIPQTD ESTLMT

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant BACE-1 protein

For Research Use Only. Not for use in diagnostic procedures.