

Human JAM-B/VE-JAM Recombinant Protein

Catalog Number: PROTP57087

Overview

Product Name	Human JAM-B/VE-JAM Recombinant Protein
Description	Human JAM-B/VE-JAM Recombinant Protein expressed in Baculovirus with hlgG-His-tag. Sequence domain: 29-238aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 50% glycerol, 1mM DTT

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P57087, 29-238aa; ADPFSAPKDQ QVVTAVEYQE AILACKTPKK TVSSRLEWKK LGRSVSFVYY QQTLQGDFKN RAEMIDFNIR
IKNVTRSDAG KYRCEVSAPS EQGQNLEEDT VTLEVLVAPA VPSCEVPSSA LSGTVVELRC QDKEGNPAPE YTWFKDGIRL LENPRLGSQS
TNSSYTMNTK TGTLQFNTVS KLDTGEYSCE ARNSVGYRRC PGKRMQVDDL NISLEPKSCD KTHTCP CPA PELLGGPSVF LFPPKPKDTL
MISRTPEVTC VVDVSHEDP EVKFNWYVDG VEVHNAKTKP REEQYNSTYR VVSVLTVLHQ DWLNGKEYKC KVS NKALPAP IEKTISKAKG
QPREPQVYTL PPSRDELTKN QVSLTCLVKG FYPDSIAVEW ESNGQPENNY KTTTPVLDSD GSFFLYSKLT VDKSRWQQGN VFSCSV MHEA
LHNHYTQKSL SLSPGKHHHH HH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human JAM-B/VE-JAM Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.