

Human PTEN Recombinant Protein

Catalog Number: PROTP60484

Overview

Product Name	Human PTEN Recombinant Protein
Description	Human PTEN Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-403aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM EDTA, 2mM DTT, 100mM NaCl, and 20% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P60484, 1-403aa; MTAIIKEIVS RNKRRYQEDG FDLDLTYIYP NIIAMGFPAE RLEGVYRNNI DDVVRFLDSK HKNHYKIYNL CAERHYDTAK FNCRVAQYPF EDHNPPQLEL IKPFCEDLDQ WLSEDDNHVA AIHCKAGKGR TGVMICAYLL HRGKFLKAQE ALDFYGEVRT RDKKGVITIPS QRRYVYYYSY LLKNHLDYRP VALLFHKMMF ETIPMFSGGT CNPQFVVCQL KVKIYSSNSG PTRREDKFMV FEFQPLPVC GDIKVEFFHK QNKMLKKDKM FHFVWNTFFI PGPEETSEKV ENGLCDQEI DSICSIERAD NDKEYLVLT TLTKNDLDKANK DKANRYFSPN FVKVLYFTKT VEEPSNPEAS SSTSVPDVS DNEPDHYRYS DTTDSDPENE PFDEDQHTQI TKV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PTEN Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.