

Human PP4 Catalytic Subunit Recombinant Protein

Catalog Number: PROTP60510

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human PP4 Catalytic Subunit Recombinant Protein
Description	Human PP4 Catalytic Subunit Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-307aa. Application(s): SDS-PAGE, Denatured. Denatured protein is most likely not the best option for functional studies. It is better suited for Western Blot (WB) or imaging assays.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P60510, 1-307aa; MGSSHHHHHH SSGLVPRGSH MGSM AEISDL DRQIEQLRRC ELIKESEVKA LCAKAREILV EESNVQRVDS
PVTVCGDIHG QFYDLKELFR VGGDVPETNY LFMGDFVDRG FYSVETFLLL LALKVRYPDRL IRLRGNHES RQITQVYGFY DECLRKYGSV
TVWRYCTEIF DYLSLSAIIID GKIFCVHGG LSPSIQTLDQI RTIDRKQEV P HDGPMCDLLW SDPEDTTGWG VSPRGAGYLF GSDVVAQFNA
ANDIDMICRA HQLVMEGYKW HFNETVLTWV SAPNYCYRCG NVAAILLEDE HLQKDFIIFE AAPQETRGIP SKKPVADYFL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human PP4 Catalytic Subunit Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.