

Human SNAP25 Recombinant Protein

Catalog Number: PROTP60880

Overview

| | |
|--------------|---|
| Product Name | Human SNAP25 Recombinant Protein |
| Description | Human SNAP25 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-206aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 50mM HEPES buffer (pH 7.4) containing 5mM DTT, 100mM NaCl |

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P60880, 1-206aa; MAEDADMRNE LEEMQRRADQ LADESLESTR RMLQLVEESK DAGIRTLVML DEQGEQLERI EEGMDQINKD MKEAEKNLTD LGKFCGLCVC PCNKLKSSDA YKKA WGNNQD GVVASQPARV VDEREQMAIS GGFIRRV TND ARENEMDENL EQVSGIIGNL RHMALDMGNE IDTQNRQIDR IMEKADSNKT RIDEANQRAT KMLGSG

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SNAP25 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.