

Human RPS3A Recombinant Protein

Catalog Number: PROTP61247

Overview

Product Name	Human RPS3A Recombinant Protein
Description	Human RPS3A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-264aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 2mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P61247, 1-264aa; MAVGKN KRLTKGGKKG AKKKVVDVDF KKDWDYDKAP AMFNIRNIGK TLVTRTQGTK IASDGLKGRV FEVSLADLQN DEVAFRKFKL ITEDVQGKNC LTNFHGMDLT RDKMCSMVKK WQTMIEAHVD VKTTDGYLLR LFCVGFTHKR NNQIRKTSYA QHQVVRQIRK KMMEIMTREV QTNDLKEVVN KLIPDSIGKD IEKACQSIYP LHDVFVRKVK MLKKPKFELG KLMELHGEYS SSGKATGDET GAKVERADGY EPPVQESV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RPS3A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.