

## Human PA28 Activator gamma Subunit Recombinant Protein

Catalog Number: PROTP61289

### Overview

Product Name	Human PA28 Activator gamma Subunit Recombinant Protein
Description	Human PA28 Activator gamma Subunit Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-254aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 40% glycerol, 200mM NaCl

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P61289, 1-254aa; MASLLKVDQE VKLKVDSFRE RITSEAEDLV ANFFPKLLE LDSFLKEPIL NIHDLTQIHS DMNLPVPDPI LLTNSHDGLD  
GPTYKKRRLD ECEEAFQGTK VVMPNGMLK SNQQLVDIIE KVKPEIRLLI EKCNTVKMWV QLLIPRIEDG NNFVSIQEE TVAELRTVES  
EAASYLDQIS RYYITRAKLV SKIAKYPHVE DYRRTVTEID EKEYISLRLLI ISELRNQYVT LHDMLKNIE KIKRPRSSNA ETLV

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PA28 Activator gamma Subunit Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.