

Human PP1 Catalytic Recombinant Protein

Catalog Number: PROTP62136

Overview

Product Name	Human PP1 Catalytic Recombinant Protein
Description	Human PP1 Catalytic Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-330aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	50mM Tris-HCl buffer (pH 8.5) containing 0.2M NaCl, 1mM DTT, 0.1mM PMSF, 1mM MnCl ₂ , 50%glycerol

Concentration

0.2mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 80% by SDS-PAGE

Amino Acid Sequence

Human, P62136, 1-330aa; MSDSEKLNLD SIIGRLLEVQ GSRPGKNVQL TENEIRGLCL KSREIFLSQP ILLELEAPLK ICGDIHGQYY DLLRLFYGG
FPPESNYLFL GDYVDRGKQS LETICLLAY KIKYPENFFL LRGNHECASI NRIYGFYDEC KRRYNIKLWK TFTDCFNCLP IAAIVDEKIF CCHGGLSPDL
QSMEQIRRIM RPTDVPDQGL LCDLLWSDPD KDVQGWGEND RGVSFTEGAE VVAKFLHKHD LDLICRAHQV VEDGYEFFAK RQLVTLFSAP
NYCGEFDNAG AMMSVDETLM CSFQILKPAD KNKGKYGQFS GLNPGGRPIT PPRNSAKAKK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PP1 Catalytic Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.