

Human RPS14 Recombinant Protein

Catalog Number: PROTP62263

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human RPS14 Recombinant Protein
Description	Human RPS14 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-151aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 40% glycerol, 100mM NaCl

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P62263, 1-151aa; MAPRK GK EKKEEQVISL GPQVAEGENV FGVCHIFASF NDTFVHVTDL SGKETICRVT GGMKVKADRD ESSPYAAMLA AQDVAQRCKE LGITALHIKL RATGGNRTKT PGPGAQSALR ALARSGMKIG RIEDVTPIPS DSTRKGGRR GRRL

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human RPS14 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.