

Human RPL30 Recombinant Protein

Catalog Number: PROTP62888

Overview

Product Name	Human RPL30 Recombinant Protein
Description	Human RPL30 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-115aa. Application(s): SDS-PAGE.
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 40% glycerol, 2mM DTT

Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

MGSSHHHHHH SGLVPRGSH MGSMVAAKKT KKSLESINSR LQLVMKSGKY VLGKQTLKM IRQGKAKLVI LANNC PALRK SEIEYYAMLA
KTVVHHYSGN NIELGTACGK YYRVCTLAII DPGDSDIIRS MPEQTGEK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RPL30 Recombinant Protein