

Mouse Vamp2 Recombinant Protein

Catalog Number: PROTP63044

Overview

| | |
|--------------|---|
| Product Name | Mouse Vamp2 Recombinant Protein |
| Description | Mouse Vamp2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-94aa. Application(s): SDS-PAGE. |
| Source | E. coli |
| Formulation | Phosphate-Buffered Saline (pH 7.4) containing 1mM EDTA, 0.1mM PMSF, 10% glycerol |

Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

MSATAA TVPPAAPAGE GGPPAPPPNL TSNRRLLQQTQ AQVDEVVDIM RVNVDPKVLER DQKLSLDDR ADALQAGASQ FETSAAKLKR
KYWWKNLK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Vamp2 Recombinant Protein