

Human RACK1 Recombinant Protein

Catalog Number: PROTP63244

Overview

Product Name	Human RACK1 Recombinant Protein
Description	Human RACK1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-317aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 50% glycerol, 0.2M NaCl, 5mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P63244, 1-317aa; MTEQMTLRGT LKGHNGWVTQ IATTPQFPDM ILSASRDKTI IMWKLTRDET NYGIPQRALR GHSHFVSDVV
ISSDGQFALS GSWDGLRLW DLTTGTTTRR FVGHTKDVLS VAFSSDNRQI VSGSRDKTIK LWNTLGVCKY TVQDESHSEW VSCVRFSPNS
SNPIIVSCGW DKLVKVVNLA NCKLKTNHIG HTGYLNTVTV SPDGSLCASG GKDGQAMLWD LNEGKHLTYL DGGDIINALC FSPNRYWLCA
ATGPSIKIWD LEGKIIVDEL KQEVISTSSK AEPQCTSLA WSADGQTLFA GYTDNLVRVW QVTIGTR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RACK1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.