

Human YBX1 Recombinant Protein

Catalog Number: PROTP67809

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human YBX1 Recombinant Protein
Description	Human YBX1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-324aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P67809, 1-324aa; MSSEAET QQPPAAPPAA PALSAADTKP GTTGSAGSG GPGGLTSAAP AGGDKKVIAT KVLGTVKWFN VRNGYGFINR
NDTKEDVVFH QTAIKKNNPR KYLRVSGDGE TVEFDVVEGE KGAEAAANTG PGGVPVQGSK YAADRNHYRR YPRRRGPPRN YQQNYQNS
GEKNEGSESA PEGQAQRRP YRRRRFPYY MRRPYGRRPQ YSNPPVQGEV MEGADNQGAG EQGRPVRQNM YRGYRPRFRR GPPRQRQPRE
DGNEEDKENQ GDETQGQPP QRRYRRNFNY RRRRPENPKP QDGKETKAAD PPAENSSAPE AEQGGA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Human YBX1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.