

Human recombinant Casein kinase 2 beta/CSNK2B protein

Catalog Number: PROTP67870-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Casein kinase 2 beta/CSNK2B protein
Description	Human recombinant Casein kinase 2 beta/CSNK2B protein expressed in E.coli with Non-Tagged. Sequence domain: 1-215aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E.coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.2M NaCl, 1mM DTT, 1mM EDTA

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P67870, 1-215aa; MSSSEEVSWI SWFCGLRGNE FFCEVDEDYI QDKFNLTGLN EQVPHYRQAL DMILDLEPDE ELEDNPNQSD
LIEQAAEMLY GLIHARYILT NRGIAQMLEK YQQGDFGYCP RVCENQPML PIGLSDIPGE AMVKLYCPKC MDVYTPKSSR HHHTDGAYFG
TGFPHMLFMV HPEYRPKRPA NQFVPRLYGF KIHPMAYQLQ LQAASNFKSP VKTIR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant Casein kinase 2 beta/CSNK2B protein](#)

For Research Use Only. Not for use in diagnostic procedures.