

Human Casein Kinase 2 alpha Recombinant Protein

Catalog Number: PROTP68400

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human Casein Kinase 2 alpha Recombinant Protein
Description	Human Casein Kinase 2 alpha Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-391aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT, 0.1mM PMSF

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P68400, 1-391aa; MGSSHHHHHH SGLVPRGSH MGSMMMSGPV PSRARVYTDV NTHRPREYWD YESHVVEWGN QDDYQLVRKL GRGKYSEVFE AINITNNEKV VVKILKPVKK KIKREIKIL ENLRGGPNII TLADIVKDPV SRTPALVFEH VNNTDFKQLY QTLTDYDIRF YMYEILKALD YCHSMGIMHR DVKPHNVIMID HEHRKRLRID WGLAEFYHPG QEYNVRVASR YFKGPELLVD YQMYDYSLDM WSLGCMLASM IFRKEPFFHG HDNYDQLVRI AKVLGTEDLY DYIDKYNIEL DPRFNDILGR HSRKRWERFV HSENQHLVSP EALDFLDKLL RYDHQSRLTA REAMEHPYFY TVVKDQARMG SSSMPGGSTP VSSANMMSGI SSVPTPSPLG PLAGSPVIAA ANPLGMPVPA AAGAQQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human Casein Kinase 2 alpha Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.