

PAFAH1B2 (NM_002572) Human Recombinant Protein

Catalog Number: PROTP68402

Overview

Product Name	PAFAH1B2 (NM_002572) Human Recombinant Protein
Description	Recombinant protein of human platelet-activating factor acetylhydrolase, isoform Ib, beta subunit 30kDa (PAFAH1B2)
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

Storage

Store at -80°C.

Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



PAFAH1B2 (NM_002572) Human Recombinant Protein