

Human RPP30 Recombinant Protein

Catalog Number: PROTP78346

Overview

| | |
|--------------|--|
| Product Name | Human RPP30 Recombinant Protein |
| Description | Human RPP30 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-268aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 5mM DTT, 200mM NaCl, 1mM EDTA |

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P78346, 1-268aa; MAVFADL DLRAGSDLKA LRGLVETAAH LGYSVVAINH IVDFKEKKQE IEKPVAVSEL FTTLPIVQGK SRPIKILTRL
TIIVSDPSHC NVLRATSSRA RLYDVVAVFP KTEKLFHIAC THLDVDLVCV TVTEKLPFYF KRPPINVAID RGLAFELVYS PAIKDSTMRR YTISSALNLM
QICKGKNVII SAAERPLEI RGPYDVANLG LLFGLSESDA KAAVSTNCRA ALLHGETRKT AFGIISTVKK PRPSEGDEDC LPASKKAKCE G

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RPP30 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.