

## Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein

Catalog Number: PROTP81117

### Overview

Product Name	Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein
Description	Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 25-420aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 1mM EDTA

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, P81117, 25-420aa; MGSSHHHHHH SGLVPRGSH MVPIDVDKTK VHNTPEVNA RIEPPDTGLY YDEYLKQVIE VLETDPHFRE KLQKADIEEI RSGRLSQELD LVSHKVRTRL DELKRQEVGR LRMLIKAKLD ALQDTGMNHH LLLKQFEHLN HQNPNTFESR DLDMLIKAAT ADLEQYDRTR HEEFKKYEMM KEHERREYLK TLSEEKREE ESKFEEMKRK HEDHPKVNHP GSKDQLKEVW EETDGLDPND FDPKTFKFLH DVNNDGFLDE QELEALFTRE LEKVYNPQNA EDDMIEMEEE RLRMREHVMS EIDNNKDRLV TLEEFRLATE KKEFLEPDSW ETLDQQQLFT EDELKEYESI IAIQENELKK RAEELQKQKE DLQRQHDHLE AQKQEYHQAV QHLEQKQLQQ GIAPSGPAGE LKFEPHT

[?/description\\_after\\_attributes?](#)

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein