

Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein

Catalog Number: PROTP81117

Overview

Product Name	Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein
Description	Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 25-420aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 1mM EDTA

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P81117, 25-420aa; MGSSHHHHHH SGLVPRGSH MVPIDVDKTK VHNTPEVNA RIEPPDTGLY YDEYLKQVIE VLETDPHFRE
KLQKADIEEI RSGRLSQELD LVSHKVRTRL DELKRQEVGR LRMLIKAKLD ALQDTGMNHH LLLKQFEHLN HQNPNTFESR DLDMLIKAAT
ADLEQYDRTR HEEFKKYEMM KEHERREYLK TLSEEKREE ESKFEEMKRK HEDHPKVNHP GSKDQLKEVW EETDGLDPND FDPKTFFKLH
DVNNDGFLDE QELEALFTRE LEKVYNPQNA EDDMIEMEEE RLRMREHVMS EIDNNKDRLV TLEEFRLATE KKEFLEPDSW ETLDQQQLFT
EDELKEYESI IAIQENELKK RAEELQKQKE DLQRQHDHLE AQQQEYHQAV QHLEQKQLQQ GIAPSGPAGE LKFEPHT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Nesfatin-1/Nucleobindin-2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.