

## Human H3F3A Recombinant Protein

Catalog Number: PROTP84243

### Overview

Product Name	Human H3F3A Recombinant Protein
Description	Human H3F3A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-136aa. Application(s): SDS-PAGE, Denatured.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.4M urea, 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P84243, 1-136aa; MGSSHHHHHH SGLVPRGSH MGSMARTKQT ARKSTGGKAP RKQLATKAAR KSAPSTGGVK KPHRYRPGTV  
ALREIRRYQK STELLIRKLP FQRLVREIAQ DFKTDLRFQS AAIGALQEAS EAYLVGLFED TNLCAIHA KR VTIMPKDIQL ARRIRGERA

🔗/description\_after\_attributes🔗

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human H3F3A Recombinant Protein