

Mouse recombinant sFRP-2/SFRP2 protein

Catalog Number: PROTP97299

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant sFRP-2/SFRP2 protein
Description	Mouse recombinant sFRP-2/SFRP2 protein expressed in Baculovirus with His-Tag. Sequence domain: 25-295aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 100mM NaCl, 20% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P97299, 25-295aa; LFLFGQP DFSYKRSNCK PIPANLQLCH GIEYQNMRLP NLLGHETMKE VLEQAGAWIP LVMKQCHPDT KKFLCSLFAP VCLDDLDETI QPCHSLCVQV KDRCAPVMSA FGFPWPDMLC CDRFPQDNDL CIPLASSDHL LPATEEAPKV CEACKTKNED DNDIMETLCK NDFALKIKVK EITYINRDTK IILETKSKTI YKLNGVSEED LKKSVLWLKD SLQCTCEEMN DINAPYLVMG QKQGGLVIT SVKRWQKGQR EFKRISRSIR KLQC

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant sFRP-2/SFRP2 protein](#)

For Research Use Only. Not for use in diagnostic procedures.