

## Mouse Crmp1 Recombinant Protein

Catalog Number: PROTP97427

### Overview

Product Name	Mouse Crmp1 Recombinant Protein
Description	Mouse Crmp1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-572aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.5) containing 40% glycerol, 0.2M NaCl, 0.1mM PMSF

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P97427, 1-572aa; MGSSHHHHHH SGLVPRGSH MGSEFMSHQG KKSIPHITSD RLLIRGGRII NDDQSFYADV YLEDGLIKQI GENLIVPGGV  
KTIEANGRMV IPGGIDVNTY LQKPSQGMTS ADDFFQGTKA ALAGGTMII DHVVPEPGSS LLTSFEKWHE AADTKSCCDY SLHVDITSWY  
DGVREELEVEL VQDKGVNSFQ VYMAYKDLYQ MSDSQLYEF TFLKGLGAVI LVHAENGLI AQEQKRILEM GITGPEGHAL SRPEELEAEA  
VFRAIAIAGR INCPVYITKV MSKSAADIIA LARKKGPLVF GEPIAASLGT DGTHYWSKNW AKAAAFVTSP PLSPDPTTPD YLTSLACGD  
LQVTGSGHCP YSTAQKAVGK DNFTLIPEGV NGIEERMTVV WDKAVATGKM DENQFVAVTS TNAAKIFNLY PRKGRIAVGS DADVVIWDPD  
KMKTTAKSH KSTVEYNIFE GMECHGSPLV VISQKIVFE DGNISVSKGM GRFIPRKPFP EHLVQRVIR SKVFGLSVSV RGMVDGPVYE  
VPATPKHAAP APSAKSSPSK HQPPPIRNLH QSNFSLGAQ IDDNPRRTG HRIVAPPGGR SNITSLG

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Crmp1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.