

## Human RBM3 Recombinant Protein

Catalog Number: PROTP98179

### Overview

|              |   |
|--------------|---|
| Product Name | Human RBM3 Recombinant Protein  |
| Description  | Human RBM3 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-157aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug  |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | E. coli   |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 50% glycerol, 5mM DTT   |

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P98179, 1-157aa; MGSSHHHHHH SSGLVPRGSH MGSMSSSEEGK LFGVGLNFNT DEQALEDHFS SFGPISEVVV VKDRETQRSR  
GFGFITFTNP EHASVAMRAM NGESLDGRQI RVDHAGKSAR GTRGGGFGAH GRGRYSRGG GDQGYGSGRY YDSRPGGYGY GYGRSRDYNG  
RNQGGYDRYS GGNYRDNYDN

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human RBM3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.