

## mycobacterium tuberculosis Antigen 85A Recombinant Protein

Catalog Number: PROTP9WQP3

### Overview

Product Name	mycobacterium tuberculosis Antigen 85A Recombinant Protein
Description	mycobacterium tuberculosis Antigen 85A Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 44-338aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	20mM Tris-HCl buffer (pH 8.0) containing, 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mycobacterium tuberculosis, P9WQP3, 44-338aa; FSRPGL PVEYLQVPSP SMGRDIKVQF QSGGANSPAL YLLDGLRAQD DFSGWDINTP  
AFEWYDQSGL SVVMPVGGQS SFYSDWYQPA CGKAGCQTYK WETFLTSELP GWLQANRHVK PTGSAVVGLS MAASSALTLA IYHPQQFVYA  
GAMSGLLDPS QAMGPTLIGL AMGDAGGYKA SDMWGPKEDP AWQRNDPLLN VGKLIANNTR VVVYCGNGKP SDLGGNLPA KFLEGFVRTS  
NIKFQDAYNA GGGHNGVDFD PDSGTHSWEY WGAQLNAMKP DLQRALGATP NTGPAPQGA

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



mycobacterium tuberculosis Antigen 85A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.