

Human MAT1A Recombinant Protein

Catalog Number: PROTQ00266

Overview

Product Name	Human MAT1A Recombinant Protein
Description	Human MAT1A Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-395aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 0.1M NaCl, 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q00266, 1-395aa; MNGPVDGLCD HSLSEGVMFM TSESVGEGHP DKICDQISDA VLDAHLKQDP NAKVACETVC KTG MVLLCGE
ITSMAMVDYQ RVVRDTIKHI GYDDSAKGFDFKTCNVLVAL EQQSPDIAQC VHLDRNEEDV GAGDQGLMFG YATDETEECM PLTIILAHKL
NARMADLRRS GLLPWL RPDS KTQVTVQYMQ DNGAVIPVRI HTIVISVQHN EDITLEEMRR ALKEQVIRAV VPAKYLDEDV VYHLQPSGRF
VIGGPQGDAG VTGRKIIVDT YGGWGAHGGG AFSGKDYTKV DRSAAYAARW VAKSLVKAGL CRRVLVQVSY AIGVAEPLSI SIFTYGTSSQK
TERELLDVWH KNFDLRPGVI VRDLDLKKPI YQKTACYGHF GRSEFPWEVP RKLVF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MAT1A Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.