

## Mouse gp130/IL-6R beta Recombinant Protein

Catalog Number: PROTQ00560

### Overview

Product Name	Mouse gp130/IL-6R beta Recombinant Protein
Description	Mouse gp130/IL-6R beta Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 23-617aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

QLLEPCGYIY PEPFVVQRGS NFTAICVLKE ACLQHYYVNA SYIVWKTNHA AVPREQVTVI NRTTSSVTFT DVVLPSVQLT CNILSFGQIE  
QNVYGVMTLS GFPPDKPTNL TCIVNEGKNM LCQWDPGRET YLETNYTLKS EWATEKFPDC QSKHGTSCMV SYMPTYVNI EVWVEAENAL  
GKVSSSEINF DPVDKVKPTP PYNLSVTNSE ELSSILKLSW VSSGLGLLD LKSDIQYRTK DASTWIQVPL EDTMSPRTSF TVQDLKPFTE YVFRIRSIKD  
SGKGYWSDWS EEASGTTYED RPSRPPSFY KTNPSHGQEY RSVRLIWKAL PLSEANGKIL DYEVIQTQSK SVSQTYYVTG TELTVNLTND  
RYVASLAARN KVGKSAAVL TIPSPHVTAA YSVVNLKAFP KDNLLWVEWT PPPKPVSKYI LEWCVLSENA PCVEDWQQED ATVNRTHLRG  
RLLESKCYQI TVTPVFATGP GGSESLKAYL KQAAPARGPT VRTKKVGKNE AVLAWDQIPV DDQNGFIRNY SISYRTSVGK EMVVHVDSSH  
TEYTLSSLSS DTLYMVRMAA YTDEGGKDG P EFTFTTPKFA QGEIELEHHH HHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse gp130/IL-6R beta Recombinant Protein