

## Mouse gp130/IL-6R beta Recombinant Protein

Catalog Number: PROTQ00560

### Overview

Product Name	Mouse gp130/IL-6R beta Recombinant Protein
Description	Mouse gp130/IL-6R beta Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 23-617aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, Q00560, 23-617aa; QLLEPCGYIY PEPVQVQRGS NFTAICVLKE ACLQHYYVNA SYIWWKTNHA AVPREQVTVI NRTTSSVTFT DVVLPSVQLT CNILSFGQIE QNVYGVMTLS GFPPDKPTNL TCIVNEGKNM LCQWDPGRET YLETNYTLKS EWATEKFPDC QSKHGTSCMV SYMPTYVNI EVWVEAENAL GKVSSSEINF DPVDKVKPTP PYNLSVTNSE ELSSILKLSW VSSGLGGLLD LKSDIQYRTK DASTWIVPL EDTMSPRTSF TVQDLKPFTE YVFRIRSIKD SGKGYWSDWS EEASGTTYED RPSRPPSFY KTNPSHGQEQ RSVRLIWKAL PLSEANGKIL DYEVIQTQSK SVSQTIVTGT TELTVNLND RYVASLAARN KVGKSAAVL TIPSPHVTAA YSVVNLKAFP KDNLLWVEWT PPPKPVSKYI LEWCVLSENA PCVEDWQQED ATVNRTHLRG RLLESKCYQI TVTPVFATGP GGSSESLKAYL KQAAPARGPT VRTKKVKGNE AVLAWDQIPV DDQNGFIRNY SISYRTSVGK EMVVHVDSSH TEYTLSSLSS DTYLMVRMAA YTDEGGKDGPEFTFTTPKFA QGEIELEHHH HHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse gp130/IL-6R beta Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.