

Mouse Apolipoprotein A-I Recombinant Protein

Catalog Number: PROTQ00623

Overview

Product Name	Mouse Apolipoprotein A-I Recombinant Protein
Description	Mouse Apolipoprotein A-I Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 25-264aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, Q00623, 25-264aa; DEPQSQW DKVKDFANVY VDAVKDSGRD YVSQFESSL GQQLNLNLE NWDTLGSTVS QLQERLGPLT
RDFWDNLEKE TDWVRQEMNK DLEEVKQKVQ PYLDEFQKKW KEDVELYRQK VAPLGAELQE SARQKLQELQ GRLSPVAEEF RDRMRTHVDS
LRTQLAPHSE QMRESLAQRL AELKSNPTLN EYHTRAKTHL KTLGEEKARPA LEDLRHSLMP MLETLKTQVQ SVIDKASETL TAQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Apolipoprotein A-I Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.