

## Human FKBP25 Recombinant Protein

Catalog Number: PROTQ00688

### Overview

Product Name	Human FKBP25 Recombinant Protein
Description	Human FKBP25 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-224aa. Application(s): SDS-PAGE, Enzyme Activity.
Size	Starting from 20ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, Q00688, 1-224aa; MAAVPQRAW TVEQLRSEQL PKKDIIFLQ EHGSDSFLAE HKLLGNIKNV AKTANKDHLV TAYNHLFETK  
RFKGTESISK VSEQVKNVKL NEDKPKETKS EETLDEGPPK YTKSVLKKGD KTNFPKKGDV VHCWYTGTLQ DGTVFDTNIQ TSAKKNNAK  
PLSFKVGVGK VIRGWDEALL TMSKGEKARL EIEPEWAYGK KGQPDAKIPP NAKLTFEVEL VDID

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human FKBP25 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.