

Human SORD Recombinant Protein

Catalog Number: PROTQ00796

Overview

Product Name	Human SORD Recombinant Protein
Description	Human SORD Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-357aa. Application(s): Enzyme Activity, SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5mM DTT, 20% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q00796, 1-357aa; MGSSHHHHHH SGLVPRGSH MAAAAPNNL SLVVHGPGL RLENYPIPEP GPNEVLLRMH SVGICGSDVH
YWEYGRIGNF IVKKPMVLGH EASGTVEKVG SSVKHLKPGD RVAIEPGAPR ENDEFCKMGR YNLSPSIFC ATPDDGNLC RFYKHNA AFC
YKLPDNTFE EGALIEPLSV GIHACRRGGV TLGHKVLVCG AGPIGMVTL VAKAMGAAQV VVTDL SATRL SKAKEIGADL VLQISKES PQ
EIARKVEGQL GCKPEVTIEC TGAEASIQAG IYATRS GGTL VLVGLGSEMT TVPLLHAAIR EVDIKGVFRY CNTWPVAISM LASKSVNVKP
LVTHRFPLEK ALEAFETFCK GLGLKIMLKC DPSDQNP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SORD Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.