

Human recombinant IL-5R alpha/IL5RA protein

Catalog Number: PROTQ01344-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant IL-5R alpha/IL5RA protein
Description	Human recombinant IL-5R alpha/IL5RA protein expressed in Baculovirus with His-Tag. Sequence domain: 21-342aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q01344, 21-342aa; ADPDLLPDEK ISLLPPVNFT IKVTGLAQVL LQWKPNPDQE QRNVNLEYQV KINAPKEDDY ETRITESKCV
TILHKGFSAS VRTILQNDHS LCLASSWASAE LHAPPGSPGT SIVNLTCTTN TTEDNYSRLR SYQVSLHCTW LVGTDAPEDT QYFLYYRYGS
WTEECQEYSK DTLGRNIACW FPRTFILSKG RDWLAVLVNG SSKHSAIRPF DQLFALHAID QINPPLNVTA EIEGTRLSIQ WEKPVSAFPI
HCFDYEVKIH NTRNGYLQIE KLMTNAFISI IDDLISKYDVQ VRAAVSSMCR EAGLWSEWSQ PIYVGNDDEHK PLREWHHHHH H

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant IL-5R alpha/IL5RA protein](#)

For Research Use Only. Not for use in diagnostic procedures.