

Human DHODH Recombinant Protein

Catalog Number: PROTQ02127

Overview

Product Name	Human DHODH Recombinant Protein
Description	Human DHODH Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 31-395aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 0.1M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q02127, 31-395aa; MGSSHHHHHH SGLVPRGSH MGSMTGDER FYAEHLMPTL QGLLDPEAH RLAVRFTSLG LLPRARFQDS
DMLEVRVLGH KFRNPVGIAA GFDKHGEAVD GLYKMGFGFV EIGSVTPKPQ EGNPRPRVFR LPEDQAVINR YGFNSHGLSV VEHLRRLRQ
KQAKLTEDGL PLGVNLGKNK TSVDAEDYA EGVRVLGPLA DYLVVNVSSP NTAGLRSLQG KAELRRLTK VLQERDGLRR VHRPAVLVKI
APDLTSQDKE DIASVVKELG IDGLIVTNTT VSRPAGLQGA LRSETGGLSG KPLRDLSTQT IREMYALTQG RVPIIGVGGV SSGQDALEKI
RAGASLVQLY TALTFWGPPV VGKVKRELEA LLKEQGGGV TDAIGADHRR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DHODH Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.