

## Human ALDH6A1 Recombinant Protein

Catalog Number: PROTQ02252

### Overview

|              |   |
|--------------|---|
| Product Name | Human ALDH6A1 Recombinant Protein   |
| Description  | Human ALDH6A1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 34-535aa. Application(s): SDS-PAGE. |
| Size         | Starting from 4ug   |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | E. coli   |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol, 1mM DTT   |

### Concentration

0.25mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Human, Q02252, 34-535aa; SSSVPTV KLFIGGKFVE SKSDKWIDIH NPATNEVIGR VPQATKAEMD AAIASCKRAF PAWADTSVLS RQQVLLRYQQ LIKENLKEIA KLITLEQGKT LADAEGDVFR GLQVVEHACS VTSLMMGETM PSITKDMPLY SYRLPLGVCA GIAPFNFPAM IPLWMFPMAM VCGNTFLMKP SERVPGATML LAKLLQDSGA PDGTLNIIHG QHEAVNFICD HPDIKAI SFV GSNKAGEYIF ERGSRHGKRV QANMGAKNHG VVMPDANKEN TLNQLVGAAF GAAGQRCSMAL STAVLVGEAK KWLPELVEHA KNLRVNAGDQ PGADLGLPIT PQAKERVCNL IDSGTKEGAS ILLDGRKIKV KGYENGNFVG PTIISNVKPN MTCYKEEIFG PVLVVLETET LDEAIQIVNN NPYGNGTAIF TTNGATARKY AHLVDVGVG VNVPIVPLP MFSFTGSRSS FRGDTNFY GK QGIQFYTLK TITSQWKEED ATLSPPAVVM PTMGR

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human ALDH6A1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.