

Mouse Angiotensin-1 receptor Recombinant Protein

Catalog Number: PROTQ02858

Overview

Product Name	Mouse Angiotensin-1 receptor Recombinant Protein
Description	Mouse Angiotensin-1 receptor Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 23-746aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, Q02858, 23-746aa; AMDLILINSL PLVSDAETSL TCIASGWHPH EPITIGRDFE ALMNQHQPDL EVTQDVTREW AKKVVWKREK
ASKINGAYFC EGRVRGQAIR IRTMKMRQQA SFLPATLMT VDRGDNVNIS FKKVLIKEED AVIYKNGSFI HSVPRHEVPD ILEVHLPDAQ
PQDAGVYSAR YIGGNLFTSA FTRLIVRRCE AQKWGPDCSR PCTTCKNNGV CHEDTGECIC PPGFMGRTCE KACEPHTFGR TCKERCSGPE
GCKSYVFLCP DPYGCSCATG WRGLQCNEAC PSGYYGPDCK LRCHCTNEEI CDRFQGCLCS QGWQGLQCEK EGRPRMTPQI EDLPDHIEVN
SGKFNPICKA SGWPLPTSEE MTLVKPDGTV LQPNDFNFTD RFSVAIFTVN RVLPPDSGVW VCSVNTVAGM VEKPFNISVK VLPEPLHAPN
VIDTGHNFAI INISSEPYFG DGPIKSKKLF YKPVNQAWKY IEVTNEIFTL NYLEPRTDYE LCVQLARPGE GGEGHPGPV RFTTASIGLP PPRGLSLLPK
SQTALNLTWQ PIFTNSEDEF YVEVERRSLQ TTSDQQNIKV PGNLTSVLLS NLVPREQYTV RARVNTKAQG EWSEELRAWT LSDILPPQPE
NIKISNITDS TAMVSWTIVD GYSSIIIR YKVQGKNEQ HIDVKIKNAT VTQYQLKGLE PETTYHVDIF AENNIGSSNP AFSHELRTLP HSPASADLGG
GKML

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Angiotensin-1 receptor Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.