

Human Aminoacylase/ACY1 Recombinant Protein

Catalog Number: PROTQ03154

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human Aminoacylase/ACY1 Recombinant Protein
Description	Human Aminoacylase/ACY1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-408aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q03154, 1-408aa; TSKGP EEEH PSVTLFRQYL RIRTVQPKPD YGA AVAFEE TARQLGLGCQ KVEVAPGYVV TVLTWPGTNP TLSSILLNSH TDVVPVFKEH WSHDPFEAFK DSEGYIYARG AQDMKCVSIQ YLEAVRRLKV EGHFRPRTIH MTFVPDEEVG GHQGMELFVQ RPEFHALRAG FALDEGIANP TDAFTVFYSE RSPWWWVRVTS TGRPGHASRF MEDTAAEKLH KVVNSILAFR EKEWQRLQSN PHLKEGSVTS VNLTKLEGGV AYNVIPATMS ASFDFRVAPD VDFKA FEEQL QSWCQAAGEG VTLEFAQKWM HPQVTP TDDS NPWWAAFSRV CKDMNLTLEP EIMPAATDNR YIRAVGV PAL GFSPMNRTPV LLHDHDERLH EAVFLRGVDI YTRLLPALAS VPALPSDS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human Aminoacylase/ACY1 Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.