

Human Glyoxalase I Recombinant Protein

Catalog Number: PROTQ04760

Overview

Product Name	Human Glyoxalase I Recombinant Protein
Description	Human Glyoxalase I Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-184aa. Application(s): Enzyme Activity, SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, Q04760, 1-184aa; ARPCIPKSFQ YSSVVCVCNA TYCDSFDPPT FPALGTFSTRY ESTRSGRRME LSMGPIQANH TGTGLLLTLQ PEQKFQKVKQ FGGAMTDAAL LNILALSPPA QNLLLKSYFS EEGIGYNIIR VPMASCDFSI RTYTYADTPD DFQLHNFSLP EEDTKLKIPL IHRALQLAQR PVSLASPWT SPTWLKTNGA VNGKGS LKGQ PGDIYHQTWA RYFVKFLDAY AEHLKQFWAV TAENEPSAGL LSGYPFQCLG FTPEHQ RDFI ARDLGPTLAN STHHNVRLLM LDDQRLLLP WAKVVLTDPE AAKYVHGI AV HWYLDLAPA KATLGETHRL FPNTMLFASE ACVGSKFWEQ SVRLGSWDRG MQYSHSIITN LLYHVVGWTD WNLALNPEGG PNWVRNFVDS PIIVDITKDT FYKQPMFYHL GHFSKFIPEG SQRVGLVASQ KNDLDAVALM HPDGSVAVVV LNRSSKDVPL TIKDPAVGFL ETISPGYSIH TYLWRRQHHH HHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Glyoxalase I Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.