

# cleavage stimulation factor (CSTF1) (NM\_001324) Human Recombinant Protein

Catalog Number: PROTQ05048

## Overview

Product Name	cleavage stimulation factor (CSTF1) (NM_001324) Human Recombinant Protein
Description	Recombinant protein of human cleavage stimulation factor, 3' pre-RNA, subunit 1, 50kDa (CSTF1), transcript variant 2
Size	20 µg
Tag	C-Myc/DDK
Form	Frozen Solution in PBS Buffer
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

## Concentration

>50 µg/mL as determined by microplate BCA method

## Storage

Store at -80°C. Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. (Ship on dry ice.)

## Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

## Amino Acid Sequence

MYRTKVGLKDRQQLYKLIISQLLYDGYISIANGLINEIKPQSVCAPSEQLLHLIKLGMENDDTAVQYAIGRSDTVAPGTGIDLEFDADVQTMSP EASEYETC  
YVTSHKGPCRVATYSRDGQLIATGSADASIKILDTERMLAKSAMPPIEVMMNETAQQNMENHPVIRTLYDHVDEVTCCLAFHPTEQILASGSRDYTLKLFDY  
SKPSAKRAFKYIQEAEMLRISIFHPSGDFILVGTQHPTLRLYDINTFQCFVSCNPQDQHTDAICSVNYNSSANMYVTGSKDGCIKLWDGVSNRCITTFEK  
AHDGAEVCSAIFSKNSKYILSSGKDSVAKLWEISTGRTLVRYTGAGLSGRQVHRTQAVFNHTEDYVLLPDERTISLCCWDSRTAERRNLLSLGHNNIVRCI  
VHSPTNPGFMTCSDDFRARFWYRRSTTD

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



cleavage stimulation factor (CSTF1) (NM\_001324) Human Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.