

Human GFPT1 Recombinant Protein

Catalog Number: PROTQ06210

Overview

Product Name	Human GFPT1 Recombinant Protein
Description	Human GFPT1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 332-699aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.15M NaCl, 50% glycerol, 2mM DTT, 2mM EDTA

Concentration

1mg/ml (determined by BRADFORD assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q06210, 332-699aa; MGSSHHHHHH SGLVPRGSH MGSQQIMKGN FSSFMQKEIF EQPESVVNTM RGRVNFDDYT VNLGGLKDHI
KEIQRCRRLI LIACGTSYHA GVATRQVLEE LTELPMVEL ASDFLDRNTP VFRDDVCFLL SQSGETADTL MGLRYCKERG ALTVGITNTV
GSSISRETDC GVHINAGPEI GVASTKAYTS QFVSLVMFAL MMCCRISMQ ERRKEIMLGL KRLPDLIKEV LSMDDDEIQKL ATELYHQKSV
LIMGRGYHYA TCLEGALKIK EITYMHSEGI LAGELKHGPL ALVDKLMPLVI MIIMRDHTYA KCQNALQQVW ARQGRPVIC DKEDTETIKN
TKRTIKVPHS VDCLQGILSV IPLQLLAFHL AVLRGYDVDF PRNLAKSVTV E

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GFPT1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.