

## Human PPP3CA Recombinant Protein

Catalog Number: PROTQ08209

### Overview

Product Name	Human PPP3CA Recombinant Protein
Description	Human PPP3CA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-511aa. Application(s): SDS-PAGE.
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5mM DTT, 1mM EDTA, 20% glycerol

### Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

MGSSHHHHHH SGLVPRGSH TGSMSPEKAI DPKLSTTDRV VKAVPFPPSH RLTAKEVFDN DGKPRVDILK AHLMKEGRLE ESVALRIITE  
GASILRQEK LLDIDAPVTV CGDIHGQFFD LMKLFEVGG S PANTRYLFLG DYVDRGYFSI ECVLYLWALK ILYPKTLFLL RGNHECRHLT  
EYFTFKQECK IKYSERVYDA CMDAFDCLPL AALMNQQFLC VHGGLSPEIN TLDDIRKLDR FKEPPAYGPM CDILWSDPLE DFGNEKTQEH  
FTHTVVRGCS YFYSYPAVCE FLQHNNLLSI LRAHEAQDAG YRMYRKSQTT GFPSLITIFS APNYLDVYNN KAAVLKYENN VMNIRQFNCS  
PHPYWLPNFM DVFTWLSLPFV GEKVTEMLVN VLNICSDDEL GSEEDGFDGA TAAARKEVIR NKIRAIGKMA RVFSVLREES ESVLTLKGLT  
PTGMPLSGVL SGGKQTLQSA IKGFSQHKI TSFEEAKGLD RINERMPPRR DAMPSDANLN SINKALTSET NGTDSNGSNS SNIQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PPP3CA Recombinant Protein