

Human Cyclophilin-40/PPiD Recombinant Protein

Catalog Number: PROTQ08752

Overview

Product Name	Human Cyclophilin-40/PPiD Recombinant Protein
Description	Human Cyclophilin-40/PPiD Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-370aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, Q08752, 1-370aa; MSHPSQAKP SNPSNPRVFF DVDIGGERVG RIVLELFADI VPKTAENFRA LCTGEGKIGH TTGKPLHFKG CPFHRIKKF MIQGGDFSNQ NGTGGESIYG EKFEDEFHY KHDREGLISM ANAGRNTNGS QFFITTVPTP HLDGKHVVFG QVIKIGVAR ILENVEVKGE KPAKLCVIAE CGELKEGDDG GIFPKDGS GD SHPDFPEDAD IDLKDVDKIL LITEDLKNIG NTFKFSQNWE MAIKKYAEVL RYVDSSKAVI ETADRKLQP IALSCVLNIG ACKLKMSNWQ GAIDSCLEAL ELDPSNTKAL YRRAQGWQGL KEYDQALADL KKAQGIAPED KAIQAELLKV KQKIKAKDK EKAVYAKMFA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Cyclophilin-40/PPID Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.